

# **Certificate of Analysis**

Laboratory Sample ID: TE40924004-019



Sep 30, 2024 | Project Packs License # 00000084ESFH12297246 2239 N Black Canyon Hwy Phoenix, AZ, 85009, US

# **Kaycha Labs**

MADH240612 Mad Honey

Classification: Hybrid

Matrix: Flower Type: Cannabis Flower

> **Production Method: Cured** Batch#: MADH240612

> > **Harvest Date:** 09/03/24

Sample Size Received: 22.02 gram

Total Amount: 7 gram

Retail Product Size: 10 gram Retail Serving Size: 10 gram

Servings: 1

Ordered: 09/24/24 Sampled: 09/24/24

Sample Collection Time: 10:30 AM

Completed: 09/27/24 Revision Date: 09/30/24

**PASSED** 

Pages 1 of 6

**SAFETY RESULTS** 







Heavy Metals **PASSED** 



Microbials **PASSED** 



**PASSED** 



Solvents **NOT TESTED** 



**NOT TESTED** 



Water Activity **NOT TESTED** 



Moisture **NOT TESTED** 



Terpenes **TESTED** 

**PASSED** 



## Cannabinoid

**Total THC** 



**Total CBD** 



**Total Cannabinoids** 



Analysis Method : SOP.T.30.500, SOP.T.30.031, SOP.T.40.031

Analytical Batch : TE005917POT Instrument Used : TE-004 "Duke Leto" (Flower) Analyzed Date : 09/24/24 19:25:42

Dilution: 400 Reagent: N/A Consumables: N/A Pipette: N/A

Analyzed by: 432, 312, 272, 87

Reviewed On: 09/26/24 12:46:29 Batch Date: 09/24/24 12:14:54

Full spectrum cannabinoid analysis utilizing High Performance Liquid Chromatography with Photo Diode Array detector (HPLC-PDA) for analysis. (Methods: SOP.T.30.500 for sample homogenization, SOP.T.30.031 for sample prep, SOP.T.40.031 for analysis on Shimadzu LC-20X0 series HPLCs). Potency results for cannabis flower products are reported on an "as received" basis. without moisture correction.

Extraction date:

09/25/24 11:33:17

This Kaycha Labs Certification shall not be reproduced, unless in its entirety, without written approval from Kaycha Labs. The results relate only to the material or product analyzed. ND=Not Detected, ppm=Parts Per Million, ppb=Parts Per Billion, RSD=Relative Standard Deviation. Limit of Detection (LOD) and Limit Of Quantitation (LOQ) are terms used to describe the smallest concentration that can be detected and reliably measured by an analytical procedure, respectively. Action Levels are State-determined thresholds based on the action limits published in Table 3.1 of 9 A.A.C. 17 and 9 A.A.C. 18. The Measurement of Uncertainty (MU) error is available from the lab upon request. The "Decision Rule" for pass/fail does not include the MU. Any calculated totals may contain rounding errors. Testing results were obtained according to requirements stated in QMS.100.010.AZ Quality Manual

Weight: 0.1949g

## **Ariel Gonzales**

Lab Director

00000024LCMD66604568 ISO 17025 Accreditation # 97164



Signature 09/27/24

Extracted by:



# Kaycha Labs

MADH240612 Mad Honey

Matrix : Flower Type: Cannabis Flower



# **Certificate of Analysis**

Project Packs

2239 N Black Canyon Hwy Phoenix, AZ, 85009, US **Telephone:** (530) 514-0500 **Email:** adam@projectpacks.co **License #:** 00000084ESFH12297246 Sample: TE40924004-019 Batch#: MADH240612

Batch#: MADH240612 Sampled: 09/24/24 Ordered: 09/24/24 Sample Size Received : 22.02 gram

Total Amount: 7 gram
Completed: 09/27/24 Expires: 09/30/25

Completed: 09/27/24 Expires: 09/30/2 Sample Method: SOP Client Method **PASSED** 

Page 2 of 6



# **Terpenes**

**TESTED** 

Terpenes	LOQ (%)	mg/g	%	Result (%)	Terpenes		LOQ (%)	mg/g	%	Result (%)
TOTAL TERPENES	0.0020	13.220	1.3220		ALPHA-CEDRENE		0.0020	ND	ND	
IMONENE	0.0020	3.792	0.3792		ALPHA-PHELLANDREN	E	0.0020	ND	ND	
BETA-CARYOPHYLLENE	0.0020	2.982	0.2982		ALPHA-TERPINENE		0.0020	ND	ND	
BETA-MYRCENE	0.0020	1.589	0.1589		ALPHA-TERPINEOL		0.0020	ND	ND	
ALPHA-HUMULENE	0.0020	1.196	0.1196		CIS-NEROLIDOL		0.0020	ND	ND	
ALPHA-PINENE	0.0020	1.080	0.1080		GAMMA-TERPINENE		0.0020	ND	ND	
INALOOL	0.0020	0.958	0.0958		GAMMA-TERPINEOL		0.0020	ND	ND	
CIMENE	0.0020	0.872	0.0872		TRANS-NEROLIDOL		0.0020	ND	ND	
BETA-PINENE	0.0020	0.751	0.0751		Analyzed by:	Weight:	Ext	raction o	date:	Extracted by:
B-CARENE	0.0020	ND	ND		334, 272, 87	0.2532g	09,	/24/24 19	9:19:39	409
ORNEOL	0.0020	ND	ND		Analysis Method : SOP.T.		0.064, SC	P.T.40.0	164	
AMPHENE	0.0020	ND	ND		Analytical Batch : TE0059 Instrument Used : TE-096		1 " TE 00	7 11 1 1 1	Farnanas	Reviewed On: 09/26/24 16:33:4 Batch Date: 09/24/24 17:26:23
AMPHOR	0.0020	ND	ND		1",TE-093 "GC - Terpene		5 1", IE-US	97 "AS -	erpenes	Batch Date : 09/24/24 17:26:23
ARYOPHYLLENE OXIDE	0.0020	ND	ND		Analyzed Date: 09/25/24					
EDROL	0.0020	ND	ND		Dilution : N/A					
UCALYPTOL	0.0020	ND	ND		Reagent: 101723.21; 07		04020-0	0000014	62. 2024	10202-1- CD22001-17215771
ENCHONE	0.0020	ND	ND		Pipette: N/A	H109203-1; 043	04030; 8	0000314	63; 2024	10202; 1; GD23001; 17315771
ENCHYL ALCOHOL	0.0020	ND	ND			med using GC-MS	which can	detect he	low single	e digit ppm concentrations. (Methods:
ERANIOL	0.0020	ND	ND		SOP.T.30.500 for sample ho	mogenization, SOF	T.30.064	for sample	e prep, an	d SOP.T.40.064 for analysis via ThermoScientifi
ERANYL ACETATE	0.0020	ND	ND							and detection carried out by ISQ 7000-series gresult is for informational purposes only and
UAIOL	0.0020	ND	ND		cannot be used to satisfy dis	pensary testing re	equirement	s in R9-17	7-317.01(A	A) or labeling requirements in R9-17-317. Nor,
OBORNEOL	0.0020	ND	ND		R9-18-310 - Q3.	ijuana establishm	ent testing	requirem	ents in R9	I-18-311(A) or labeling requirements in
SOPULEGOL	0.0020	ND	ND							
IENTHOL	0.0020	ND	ND							
EROL	0.0020	ND	ND							
ULEGONE	0.0020	ND	ND							
ABINENE	0.0020	ND	ND							
ABINENE HYDRATE	0.0020	ND	ND							
ERPINOLENE	0.0020	ND	ND							
ALENCENE	0.0020	ND	ND							
ALPHA-BISABOLOL	0.0020	NID	ND							

This Kaycha Labs Certification shall not be reproduced, unless in its entirety, without written approval from Kaycha Labs. The results relate only to the material or product analyzed. ND=Not Detected, ppm=Parts Per Million, ppb=Parts Per Billion, RSD=Relative Standard Deviation. Limit of Detection (LOD) and Limit Of Quantitation (LOQ) are terms used to describe the smallest concentration that can be detected and reliably measured by an analytical procedure, respectively. Action Levels are State-determined thresholds based on the action limits published in Table 3.1 of 9 A.A.C. 17 and 9 A.A.C. 18. The Measurement of Uncertainty (MU) error is available from the lab upon request. The "Decision Rule" for pass/fail does not include the MU. Any calculated totals may contain rounding errors. Testing results were obtained according to requirements stated in QMS.100.010.AZ Quality Manual.

## **Ariel Gonzales**

Lab Director

State License # 00000024LCMD66604568 ISO 17025 Accreditation # 97164 ait Dong

Signature 09/27/24



# **Kaycha Labs**

MADH240612 Mad Honey Matrix: Flower

Type: Cannabis Flower

# **Certificate of Analysis**

**PASSED** 

2239 N Black Canyon Hwy Phoenix, AZ, 85009, US Telephone: (530) 514-0500 Email: adam@proiectpacks.co **License # :** 00000084ESFH12297246 Sample: TE40924004-019 Batch#: MADH240612 Sampled: 09/24/24 Ordered: 09/24/24

Sample Size Received: 22.02 gram Total Amount: 7 gram
Completed: 09/27/24 Expires: 09/30/25 Sample Method: SOP Client Method

Page 3 of 6



# **Pesticides**

PAS	SE	
-----	----	--

Pesticide	LOQ	Units		rel Pass/Fail	Res
AVERMECTINS (ABAMECTIN B1A)	0.2500	ppm	0.5	PASS	ND
ACEPHATE	0.2000	ppm	0.4	PASS	ND
ACETAMIPRID	0.1000	ppm	0.2	PASS	ND
ALDICARB	0.2000	ppm	0.4	PASS	ND
AZOXYSTROBIN	0.1000	ppm	0.2	PASS	ND
BIFENAZATE	0.1000	ppm	0.2	PASS	ND
BIFENTHRIN	0.1000	ppm	0.2	PASS	ND
BOSCALID	0.2000	ppm	0.4	PASS	ND
CARBARYL	0.1000	ppm	0.2	PASS	ND
CARBOFURAN	0.1000	ppm	0.2	PASS	ND
CHLORANTRANILIPROLE	0.1000	ppm	0.2	PASS	ND
CHLORPYRIFOS	0.1000	ppm	0.2	PASS	ND
CLOFENTEZINE	0.1000	ppm	0.2	PASS	ND
CYPERMETHRIN	0.5000	ppm	1	PASS	ND
DIAZINON	0.1000	ppm	0.2	PASS	ND
DAMINOZIDE	0.5000	ppm	1	PASS	ND
DICHLORVOS (DDVP)	0.0500	ppm	0.1	PASS	ND
DIMETHOATE	0.1000	ppm	0.2	PASS	ND
ETHOPROPHOS	0.1000	ppm	0.2	PASS	ND
ETOFENPROX	0.2000	ppm	0.4	PASS	ND
ETOXAZOLE	0.1000	ppm	0.2	PASS	ND
FENOXYCARB	0.1000	ppm	0.2	PASS	ND
FENPYROXIMATE	0.2000	ppm	0.4	PASS	ND
FIPRONIL	0.2000	ppm	0.4	PASS	ND
FLONICAMID	0.5000	ppm	1	PASS	ND
FLUDIOXONIL	0.2000	ppm	0.4	PASS	ND
HEXYTHIAZOX	0.5000	ppm	1	PASS	ND
IMAZALIL	0.1000	ppm	0.2	PASS	ND
IMIDACLOPRID	0.2000	ppm	0.4	PASS	ND
KRESOXIM-METHYL	0.2000	ppm	0.4	PASS	ND
MALATHION	0.1000	ppm	0.2	PASS	ND
METALAXYL	0.1000	ppm	0.2	PASS	ND
METHIOCARB	0.1000	ppm	0.2	PASS	ND
METHOMYL	0.2000	ppm	0.4	PASS	ND
MYCLOBUTANIL	0.1000	ppm	0.2	PASS	ND
NALED	0.2500	ppm	0.5	PASS	ND
OXAMYL	0.5000	ppm	1	PASS	ND
PACLOBUTRAZOL	0.2000	ppm	0.4	PASS	ND
TOTAL PERMETHRINS	0.1000	ppm	0.2	PASS	ND
PHOSMET	0.1000	ppm	0.2	PASS	ND
PIPERONYL BUTOXIDE	1.0000	ppm	2	PASS	ND
PRALLETHRIN	0.1000	ppm	0.2	PASS	ND
PROPICONAZOLE	0.2000	ppm	0.4	PASS	ND
PROPOXUR	0.1000	ppm	0.2	PASS	ND
TOTAL PYRETHRINS	0.5000	ppm	1	PASS	ND
PYRIDABEN	0.1000	ppm	0.2	PASS	ND

Pesticide		LOQ	Units	Action Level	Pass/Fail	Result	
TOTAL SPINOSAD		0.1000	ppm	0.2	PASS	ND	
SPIROMESIFEN		0.1000	ppm	0.2	PASS	ND	
SPIROTETRAMAT		0.1000	ppm	0.2	PASS	ND	
SPIROXAMINE		0.2000	ppm	0.4	PASS	ND	
TEBUCONAZOLE		0.2000	ppm	0.4	PASS	ND	
THIACLOPRID		0.1000	ppm	0.2	PASS	ND	
THIAMETHOXAM		0.1000	ppm	0.2	PASS	ND	
TRIFLOXYSTROBIN		0.1000	ppm	0.2	PASS	ND	
CHLORFENAPYR *		0.3000	ppm	1	PASS	ND	
CYFLUTHRIN *		0.5000	ppm	1	PASS	ND	
Analyzed by: 152, 272, 87	Weight: 0.4991g	Extraction of 09/25/24 12			Extracted 410	l by:	
Analytical Batch :TE005919 Instrument Used :TE-117 "U Analyzed Date :09/25/24 19	JHPLC - Pest/Myco 1",TE-26	2 "MS/MS - Pest/Myc	0 2"	Reviewed On: 09/27/24 16:29:50 Batch Date: 09/24/24 12:18:47			
Dilution: 25 Reagent: 091324.R12; 090: Consumables: 947.155; 80 Pipette: TE-060 SN:20C354	00038072; 111423CH01; 2: 57 (20-200uL); TE-108 SN:2	20318-306-D; 10086 20B18337 (100-1000	45998; GD2: uL)	3001; 425240JF			
Pesticide screening is carried homogenization, SOP.T.30.104							
Analyzed by: 152, 272, 87	Weight: 0.4991g	Extraction 6 09/25/24 12			Extracted 410	l by:	
Analysis Method: SOP.T.30. Analytical Batch: TE005952 Instrument Used: TE-117 "I Analyzed Date: 09/27/24 12	VOL NS/MS Pest/Myco 1",TE-262		2		n:09/27/24 16:09/27/24 12:1		
Dilution: 25							

Dilution: 25
Reagent: 091324.R12; 090524.R14; 091324.R13; 073024.R30; 091924.R02; 091824.R01; 091324.R31; 091924.R03; 041823.06
Consumables: 947.155; 8000038072; 111423CH01; 220318-306-D; 1008645998; GD23001; 425240]F
Pipette: 1TE-060 SN:20C33457 (20-200L); TE-108 SN:20B18337 (100-1000L)
Supplemental pesticide screening using GC-MS/R05 to quantitatively screen for Chiofrenapyr, Cyfluthrin, Cypermethrin, and Diazinor; as well as the qualitative confirmation of Dichlorvos, Permethrins, Piperonyl Butoxide, Prailethrin, Propiconazole, Pyrethrins, and Tebuconazole which are all quantitatively screened using LC-MS/RMS. (Methods: SoPT-30.500 for sample homogenetation, SOPT-31.014.AZ for sample and SOPT-40.14AZ for analysis using a ThermoScietific 1310-series GC equipped with a TriPlus RSH autosampler and detected on a TSQ 9000-series mass spectrometer).

This Kaycha Labs Certification shall not be reproduced, unless in its entirety, without written approval from Kaycha Labs. The results relate only to the material or product analyzed. ND=Not Detected, ppm=Parts Per Million, ppb=Parts Per Billion, RSD=Relative Standard Deviation. Limit of Detection (LOD) and Limit Of Quantitation (LOQ) are terms used to describe the smallest concentration that can be detected and reliably measured by an analytical procedure, respectively. Action Levels are State-determined thresholds based on the action limits published in Table 3.1 of 9 A.A.C. 17 and 9 A.A.C. 18. The Measurement of Uncertainty (MU) error is available from the lab upon request. The "Decision Rule" for pass/fail does not include the MU. Any calculated totals may contain rounding errors. Testing results were obtained according to requirements stated in QMS.100.010.AZ Quality Manual

## **Ariel Gonzales**

Lab Director

00000024LCMD66604568 ISO 17025 Accreditation # 97164

Signature 09/27/24



# Kaycha Labs

MADH240612 Mad Honey Matrix: Flower



Type: Cannabis Flower

# **Certificate of Analysis**

PASSED

**Project Packs** 

2239 N Black Canyon Hwy Phoenix, AZ, 85009, US Telephone: (530) 514-0500 Email: adam@proiectpacks.co **License # :** 00000084ESFH12297246 Sample: TE40924004-019

Batch#: MADH240612 Sampled: 09/24/24 Ordered: 09/24/24

Sample Size Received: 22.02 gram Total Amount: 7 gram
Completed: 09/27/24 Expires: 09/30/25

Sample Method: SOP Client Method

Page 4 of 6



# Microbial



OCHRATOXIN A

# **Mycotoxins**



PASS

20

ND

Analyte		LOQ	Units	Result	Pass / Fail	Action Level
SALMONELLA S	PP	0.0000		Not Present in 1g	PASS	
ASPERGILLUS F	LAVUS	0.0000		Not Present in 1g	PASS	
ASPERGILLUS F	UMIGATUS	0.0000		Not Present in 1g	PASS	
ASPERGILLUS N	IIGER	0.0000		Not Present in 1g	PASS	
ASPERGILLUS T	ERREUS	0.0000		Not Present in 1g	PASS	
ESCHERICHIA C	OLI REC	10.0000	CFU/g	<10	PASS	100
Analyzed by: 87, 39, 272	Weight: 1.0227a	Extraction 09/25/2	on date: 4 14:51:4		Extracted 331	by:
07, 33, 272	1.02279	09/23/2	4 14.31.4	40 3	))I	

Analysis Method: SOP.T.40.056B, SOP.T.40.058.FL, SOP.T.40.208, SOP.T.40.209.AZ Analytical Batch: TE005914MIC Reviewed On: 09/26/24 11:51:03

Instrument Used : TE-234 "bioMerieux GENE-UP" **Batch Date :** 09/24/24 12:03:17

Analyzed Date : N/A

Dilution: 10

081324.13; 081324.20 Consumables: N/A Pipette: N/A

Analyte	LOQ Units	Result	Pass / Fail	Action Level
TOTAL AFLATOXINS	4.8510 ppb	ND	PASS	20
AFLATOXIN B1	4.8510 ppb	ND	PASS	20
AFLATOXIN B2	5.9400 ppb	ND	PASS	20
AFLATOXIN G1	6.2700 ppb	ND	PASS	20
AFLATOXIN G2	10.7250 ppb	ND	PASS	20

12.0000 ppb Extracted by: Weight: 0.4991g Extraction date: 09/25/24 12:41:29 Analysis Method : SOP.T.30.500, SOP.T.30.104.AZ, SOP.T.40.104.AZ

Analytical Batch : TE005951MYC Instrument Used : N/A

**Analyzed Date :** 09/27/24 12:13:02

Reviewed On: 09/27/24 16:33:07 Batch Date: 09/27/24 12:11:54

Dilution: 25 Reagent: 091624.R20; 081224.20; 081324.01; 081324.47; 081324.50; 081324.55; 081324.66; Reagent: 091324.R12; 090524.R14; 091324.R13; 073024.R30; 091924.R02; 091824.R01;

 $091324.R31; 091924.R03; 041823.06 \\ \textbf{Consumables}: 947.155; 8000038072; 111423CH01; 220318-306-D; 1008645998; GD23001; 320318-306-D; 320318-D; 320518-D; 320518-D; 320518-D; 320518$ 

425240IF

Pipette: TE-060 SN:20C35457 (20-200uL); TE-108 SN:20B18337 (100-1000uL)

Aflatoxins B1, B2, G1, G2, and Ochratoxin A analysis using LC-MS/MS. (Methods: SOP.T.30.500 for sample homogenization, SOP.T.30.104.AZ for sample prep, and SOP.T.40.104.AZ for analysis on ThermoScientific Altis TSQ with Vanquish UHPLC). Total Aflatoxins (sum of Aflotoxins B1, B2, G1, G2) must be <20 $\mu$ g/kg. Ochratoxin must be <20µg/kg



# **Heavy Metals**

# **PASSED**

Metal		LOQ Units	Result	Pass / Fail	Action Level
ARSENIC		0.2000 ppm	ND	PASS	0.4
CADMIUM		0.2000 ppm	ND	PASS	0.4
LEAD		0.5000 ppm	ND	PASS	1
MERCURY		0.6000 ppm	ND	PASS	0.2
Analyzed by:	Weight:	Extraction date:		Extracted	l by:
398, 39, 272, 87	0.191g	09/25/24 15:13:21		398	

Analysis Method: SOP.T.30.500, SOP.T.30.084.AZ, SOP.T.40.084.AZ

Analytical Batch : TE005930HEA

**Reviewed On:** 09/26/24 09:52:40

Instrument Used: TE-051 "Metals Hood", TE-141 "Wolfgang", TE-153Batch Date: 09/24/24

"Bill",TE-154 "Bill's PC",TE-157 "Bill Pump",TE-156 "Bill Chiller",TE-155 "Bill AS",TE-218 "Bill Monitor",TE-219 "Bill Monitor",TE-260 "Ludwig"

Analyzed Date : N/A

Dilution: 50

Reagent: 101723.14; 092324.R01; 091624.R19; 032724.07; 081624.01; 100121.01

Consumables: 111423CH01; 210705-306-D; 210725-598-D Pipette: TE-063 SN:20C50490 (20-200uL); TE-110 SN:20B18338 (100-1000uL)

Heavy Metals screening is performed using ICP-MS (Inductively Coupled Plasma – Mass Spectrometer) which can screen down to below single digit ppb concentrations for regulated heavy metals. (Methods: SOP.T.30.500 for sample homogenization, SOP.T.30.084.AZ for sample prep by microwave digestion, and SOP.T.40.084.AZ for analysis by ThermoScientific ICAP RQ ICP-MS).

This Kaycha Labs Certification shall not be reproduced, unless in its entirety, without written approval from Kaycha Labs. The results relate only to the material or product analyzed. ND=Not Detected, ppm=Parts Per Million, ppb=Parts Per Billion, RSD=Relative Standard Deviation. Limit of Detection (LOD) and Limit Of Quantitation (LOQ) ppulparis Per Billion, RSDE-Relative Standard Deviation. Limit of Detection (LOD) and Limit of Quantitation (LOQ) are terms used to describe the smallest concentration that can be detected and reliably measured by an analytical procedure, respectively. Action Levels are State-determined thresholds based on the action limits published in Table 3.1 of 9 A.A.C. 17 and 9 A.A.C. 18. The Measurement of Uncertainty (MU) error is available from the lab upon request. The "Decision Rule" for pass/fail does not include the MU. Any calculated totals may contain rounding errors. Testing results were obtained according to requirements stated in QMS.100.010.AZ Quality Manual

## **Ariel Gonzales**

Lab Director

00000024LCMD66604568 ISO 17025 Accreditation # 97164

Signature 09/27/24



# Kaycha Labs

MADH240612 Mad Honey Matrix : Flower



Type: Cannabis Flower

# **PASSED**

# Certificate of Analysis

Project Packs

2239 N Black Canyon Hwy Phoenix, AZ, 85009, US **Telephone:** (530) 514-0500 **Email:** adam@projectpacks.co **License #:** 00000084ESFH12297246 Sample: TE40924004-03 Batch#: MADH240612 Sampled: 09/24/24 Ordered: 09/24/24

Sample Size Received: 22.02 gram Total Amount: 7 gram Completed: 09/27/24 Expires: 09/30/25 Sample Method: SOP Client Method

Page 5 of 6

# **COMMENTS**

\* Confident Cannabis sample ID: 2409KLAZ0647.2680



\* Cannabinoid

TE40924004-019POT

**1 -** M3:D9-THC

This Kaycha Labs Certification shall not be reproduced, unless in its entirety, without written approval from Kaycha Labs. The results relate only to the material or product analyzed. ND=Not Detected, ppm=Parts Per Million, ppb=Parts Per Billion, RSD=Relative Standard Deviation. Limit of Detection (LOD) and Limit Of Quantitation (LOQ) are terms used to describe the smallest concentration that can be detected and reliably measured by an analytical procedure, respectively. Action Levels are State-determined thresholds based on the action limits published in Table 3.1 of 9 A.A.C. 17 and 9 A.A.C. 18. The Measurement of Uncertainty (MU) error is available from the lab upon request. The "Decision Rule" for pass/fail does not include the MU. Any calculated totals may contain rounding errors. Testing results were obtained according to requirements stated in QMS.100.010.AZ Quality Manual.

**Ariel Gonzales** 

Lab Director

State License # 00000024LCMD66604568 ISO 17025 Accreditation # 97164 atil Jongh

Signature 09/27/24



# Kaycha Labs

MADH240612 Mad Honey Matrix : Flower



Type: Cannabis Flower

# **PASSED**

# **Certificate of Analysis**

**Project Packs** 

2239 N Black Canyon Hwy Phoenix, AZ, 85009, US **Telephone:** (530) 514-0500 **Email:** adam@projectpacks.co **License #:** 00000084ESFH12297246 Sample: TE40924004-0: Batch#: MADH240612 Sampled: 09/24/24 Ordered: 09/24/24

Sample Size Received: 22.02 gram
Total Amount: 7 gram
Completed: 09/27/24 Expires: 09/30/25
Sample Method: SOP Client Method

Page 6 of 6

# **COMMENTS**

\* Confident Cannabis sample ID: 2409KLAZ0647.2680



This Kaycha Labs Certification shall not be reproduced, unless in its entirety, without written approval from Kaycha Labs. The results relate only to the material or product analyzed. ND=Not Detected, ppm=Parts Per Million, ppb=Parts Per Billion, RSD=Relative Standard Deviation. Limit of Detection (LOD) and Limit Of Quantitation (LOQ) are terms used to describe the smallest concentration that can be detected and reliably measured by an analytical procedure, respectively. Action Levels are State-determined thresholds based on the action limits published in Table 3.1 of 9 A.A.C. 17 and 9 A.A.C. 18. The Measurement of Uncertainty (MU) error is available from the lab upon request. The "Decision Rule" for pass/fail does not include the MU. Any calculated totals may contain rounding errors. Testing results were obtained according to requirements stated in QMS.100.010.AZ Quality Manual.

## **Ariel Gonzales**

Lab Director

State License # 00000024LCMD66604568 ISO 17025 Accreditation # 97164 atil Jensh

Signature 09/27/24